



### Enhanced performance by compression

Activate your muscles before you do sports  
 Increase your performance while you do sports  
 Help you to regenerate more quickly after doing sports  
 Stabilize ligaments and joints  
 Ergonomically right/left design  
 Hand linked seam over toes for extra comfort  
 Padded sole  
 Toe and heel cushion  
 Measuring point is the widest area of the calf  
 Calf circumference in cm: I (25-31), II (32-38), III (39-44), IV (45-50)

|                               |                                           |
|-------------------------------|-------------------------------------------|
| <b>Fabric:</b>                | Outer fabric: 85% polyamide, 15% elastane |
| <b>Country of origin:</b>     | Volksrepublik China                       |
| <b>Customs tariff number:</b> | 61159610                                  |

### Care instructions:



### Available colours

|                                                                        | I      | II     | III    | IV     |
|------------------------------------------------------------------------|--------|--------|--------|--------|
| <b>Weight in g</b>                                                     | 63 g   | 66 g   | 77 g   | 88 g   |
| <b>VPE</b><br>(Pcs. per inner packaging<br>/ pcs. per outer packaging) | 10/100 | 10/100 | 10/100 | 10/100 |

| Measurements in cm                 | I        | II       | III      | IV       |
|------------------------------------|----------|----------|----------|----------|
| <b>Size</b>                        | 38,00    | 41,00    | 44,00    | 47,00    |
| <b>Size UK</b>                     | 5,50     | 7,50     | 10,00    | 12,00    |
| <b>Circumference of calf up to</b> | 31,00 cm | 38,00 cm | 44,00 cm | 50,00 cm |

### Available colours

■ black (blackC)

■ white (white)

### Features



### OEKO-TEX® Standard 100

OEKO-Tex® CONFIDENCE IN TEXTILES STANDARD 100 15.0.70467 HOHENSTEIN HTTI Tested for harmful substances. [www.oeko-tex.com/standard100](http://www.oeko-tex.com/standard100)

